# Coracto™

## RAPAMYCIN ELUTING STENT

EN	ENGLISH.....	4
TR	TÜRKÇE.....	7
FR	FRANÇAIS.....	10
DE	DEUTSCH.....	14
NL	NEDERLANDS.....	17
ES	ESPAÑOL.....	21
IT	ITALIANO.....	25
PT	PORTUGUÊS.....	28
HU	MAGYAR.....	31
HR	HRVATSKI.....	35
BS	BOSANSKI.....	38
CS	ČEŠTINA.....	41
DA	DANISH.....	45
SR	SRPSKI.....	48
SV	SWENSKA.....	51
RU	РУССКИЙ.....	54
PL	POLSKI.....	58
RO	ROMÂNĂ.....	61
MK	МАКЕДОНСКИ.....	65
BG	БЪЛГАРСКИ.....	68
EL	ΕΛΛΗΝΙΚΑ.....	72
SK	SLOVENČINA.....	76
SLO	SLOVENŠČINA.....	79
UA	УКРАЇНСЬКА.....	82
N	NORSK.....	86
LT	LIETUVA.....	89
LV	LATVIEŠU.....	92
FIN	SUOMI.....	96
EST	EESTI.....	99
AR	العربية.....	104
FA	فارس.....	106

**Rx only**

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## ENGLISH

### CORACTO™ RAPAMYCIN-ELUTING CORONARY STENT DELIVERY SYSTEM

Carefully read all instructions prior to use. Observe all warnings and precautions noted in these instructions. Failure to do so may result in complications.

#### DEVICE DESCRIPTION

##### 1. The Coracto™ Rapamycin-Eluting Coronary Stent Delivery System includes:

The CORACTO™ coronary stent system is an EtO-sterilized, single-use product made up of a balloon-expandable intracoronary stent pre-mounted on PTCA balloon rapid exchange (MonoRail) delivery system. The stent is manufactured from 316LVM medical grade stainless steel. The stents are available in multiple lengths including 9 mm, 13 mm, 17 mm, 21 mm, 24 mm, 28 mm, 32 mm and can be expanded to different diameters including 2.5 mm, 2.75 mm, 3.0 mm, 3.5 mm, and 4.0 mm. The maximum guide wire diameter appropriate for use with the system is 0.014". The stent coating is composed of the active ingredient rapamycin and the inactive ingredient PLGA, a copolymer of poly (D, L-lactic-co-glycolic-acid in a ratio of 50:50. The drug load is 1.7 µg rapamycin/mm<sup>2</sup> total stent surface.

Stent Inner Diameter (mm)	Stent Length (mm)	Minimum Guiding Catheter Compatibility ID	Nominal Pressure (atm)	Rated Burst Pressure (atm)
2.50	9-13-17-21-24-28	5F (1.4 mm / 0.056")	8	16
2.75	9-13-17-21-24-28	5F (1.4 mm / 0.056")	8	16
3.00	9-13-17-21-24-28-32	5F (1.4 mm / 0.056")	8	16
3.50	9-13-17-21-24-28-32	5F (1.4 mm / 0.056")	8	16
4.00	9-13-17-21-24-28-32	5F (1.4 mm / 0.056")	8	16

##### 2. Drug Component Description

Rapamycin is a macrocyclic lactone product of the soil bacteria *Streptomyces hygroscopicus*. Its molecular formula is C<sub>48</sub>H<sub>79</sub>NO<sub>13</sub> and its molecular weight is 914.2g/mol. The chemical name for rapamycin is (3S, 6R, 7E, 9R, 10R, 12R, 14S, 15E, 17E, 19E, 21S, 23S, 26R, 27R,34aS)9,10,12,13,14,21,22,23,24,25,26,27,32,33,34,34a-hexadecahydro-9,27-dihydroxy-3-[(1R)-2-[(1S,3R,4R)-4-hydroxy-3-methoxycyclohexyl]-1-methylethyl]-10,21-dimethoxy-6,8,12,14,20,26-hexamethyl-23,27-epoxy-3H-pyrido[2,1-c][1,4]oxazacycloheptanone-1,5,11,28,29 (4H,6H,31H)-pentone.

#### HOW PRODUCT IS SUPPLIED

The Coracto™ Rapamycin-Eluting Coronary Stent Delivery System supplied with single Instructions for Use (IFU) in each packaging box.

**Sterile:** This device is a sterilized with Ethylene Oxide. The outer foil pouch is only intended to provide protection from light and humidity (a small paper packet containing drying agent (silica gel) is included for the absorption of residual humidity). The contents of the sterile packaging placed inside the foil pouch are sterile. Dispose of the silica gel sachet together with the packaging.

**Single use:** This device is for single use only. Do not re-sterilize. Do not use if the package is opened or damaged.

**Storage:** Store in a dry, dark, cool place. Protect from light. Do not remove from the carton until ready to use. Store at 25°C; excursions permitted to 10-30°C.

#### INDICATIONS

Patients who have symptomatic ischemic heart disease associated with the stenotic or restenotic lesions (≤30 mm) in the coronary arteries (with a reference vessel diameter of 2.5 mm to 4.0 mm) and who are eligible for PTCA (Percutaneous Transluminal Coronary Angioplasty) procedures.

#### CONTRAINDICATIONS

- Patients in whom antiplatelet and/or anticoagulation therapy is contraindicated.

- Patients judged to have a lesion which prevents complete inflation of an angioplasty balloon or proper placement of the stent or delivery system.
- Patients with hypersensitivity or contraindication to sirolimus (rapamycin) or structurally-related compounds zotarolimus, everolimus, tacrolimus, or similar drugs, or any analogue or derivative
- Patients with known allergy to metal or contrast media.

#### WARNINGS

- The product should only be used by experienced physicians who have been adequately trained for percutaneous coronary intervention techniques.
- The procedure should only be performed at hospitals where emergency coronary artery bypass graft surgery can be quickly performed in the event of potentially injurious or life-threatening complications.
- Judicious selection of patients is necessary in case of thrombosis, vascular complications and / or bleeding events risks.
- Before using the product, read all Instructions for use thoroughly prior to the procedure.
- Do not use the product where the integrity of the sterile packaging has been compromised or if the device appears damaged or kinked.
- Do not use the product after the use-by date.
- Do not re-use, re-process or re-sterilize the product. This will affect the product performance. The product may not be able to withstand the stresses of re-use. Cleaning or re-processing of the catheter components exposed to biological materials may result in adverse device effect when used.
- Appropriate anticoagulant and/or antiplatelet therapy should be applied to the patients. In case of inappropriate medical therapy, stent thrombosis may occur.
- When implanting overlapping stents, use the stents of the same material in order to prevent the possibility of metal corrosion.
- Persons allergic to stainless steel may suffer an allergic reaction to this implant.
- Should any resistance be felt at any time during either lesion access or removal of delivery system post-stent implantation, the entire system should be removed as a single unit.
- The CORACTO™ Stent reaches its nominal dilated diameter when the delivery system balloon is inflated at a pressure of 8 bars. The delivery system balloon should not be inflated at pressures so high as to dilate the stent by more than the intended vessel lumen diameter or at pressures that exceed the pre-determined maximum pressure.
- After use, this product might be a potential biohazard. Handling and disposal of the products should be in accordance with accepted medical practice and applicable local laws and regulations. Alvimedica is not responsible for hazards and/or other problems which may result from incorrect handling and/or disposal of the products.

#### PRECAUTIONS

- Do not remove the stent from the delivery system. This stent cannot be removed and placed on another balloon catheter for deployment.
- In case of any compromise to the legibility of the label that affects the proper identification of the product, it is recommended not to use the device.
- Do not expose to light for excessive amount of time.
- Keep the product in a cool, dry and dark place.
- Do not exceed the Rated Burst Pressure (RBP) as indicated on the product label when inflating the balloon. This might lead to serious complications such as dissection, perforation, damage to intima, thrombosis, myocardial ischemia, additional surgical intervention, myocardial infarction, arrhythmia, stent embolism, total occlusion of coronary artery and death.
- Do not expand the stent if it is not properly positioned in the vessel.
- Incomplete expansion of the stent may cause procedural complications resulting in patient injury.

- Do not remove the stent from its delivery system as removal may damage the stent and/or lead to stent embolization. The coronary stent delivery system is used as a single unit.
- Never advance or retract the device against unknown resistance, as this could cause tissue trauma or device damage. If unusual resistance is felt at any time during lesion access before stent implementation, the stent system and guiding catheter should be removed as a single unit.
- Ensure the balloon is fully deflated before withdrawing the balloon catheter from the treatment area. The whole procedure should be performed under high quality fluoroscopy imaging. If stent location cannot be verified, do not expand the stent.
- During the procedure appropriate anticoagulant and coronary vasodilator therapy must be provided to the patient as needed. Anticoagulant therapy should be given according to the current medical guidelines and continued for a period of time to be determined by the physician after the procedure.
- Placement of the stent has the potential to compromise side branch patency.
- When treating multiple lesions, the distal lesion should be initially stented, followed by stenting of the proximal lesion. This technique alleviates the need to cross the proximal stent in placement of the distal stent and reduces the changes for dislodgement.
- Do not use air or any gaseous medium to inflate the balloon.
- Do not use the contrast materials Ethiodol or Lipiodol.
- Do not expose the delivery system to organic solvents (i.e. alcohol, etc.).
- It is recommended to remove the peel-off label for traceability purposes before disposal of the packaging.
- Care must be exercised when crossing a newly deployed stent with an intravascular ultrasound (IVUS), a coronary guidewire, or a balloon catheter to avoid disrupting the stent placement, opposition and/or geometry.
- After use, this product might be a potential biohazard. Handling and disposal of the products in accordance with accepted medical practices and applicable local laws and regulations must be done. Handling and disposal of the products is not the responsibility of Alvimedica.

#### MRI SAFETY

Non-clinical testing has demonstrated that thStent is MR Conditional. Patients can be scanned safely under the following conditions:

- Static magnetic field of 1.5-Tesla (1.5T) or 3.0-Tesla (3.0T).
- Maximum spatial gradient field less than or equal to 10T/m.
- Normal operating mode: Maximum whole body specific absorption rate (SAR) of :
  - 2.0 W/kg for 15 minutes of scanning at normal operating mode at 1.5T.
  - 2.0 W/kg for 15 minutes of scanning at normal operating mode at 3.0T.

#### 3.0 T RF Heating

In non-clinical testing with body coil excitation, the Stent produced a maximal differential temperature rise of 1.3°C when exposed to a maximum specific absorption rate (SAR) of 3.4 W/kg for 15 minutes of scanning in a 3.0T MR system (Siemens Trio, SYNGO MR A30 4VA30A software, Munich, Germany). Scaling of the SAR and observed heating indicates that SAR of 2.0 W/kg would be expected to yield a localized temperature rise of less than 1.0 °C.

#### 1.5 T RF Heating

In non-clinical testing with body coil excitation, the Stent produced a differential temperature rise of 1.1 °C when exposed to a maximum specific absorption rate (SAR) of 1.4 W/kg for 15 minutes of scanning in a 1.5T MR system (Siemens Espree, SYNGO MR B17 software, Munich, Germany). Scaling of the SAR and observed heating indicates that SAR of 2.0 W/kg would be expected to yield a localized temperature rise of less than 2.0 °C.

**Caution:** The RF heating behavior does not scale with static field strength. Devices which do not exhibit detectable heating at one field strength may exhibit high values of localized heating at another field strength.

#### MR Artifact

MR image quality may be compromised if the area of interest is the same or relatively close to the position of the device, and it may be necessary to optimize the MR imaging parameters. The shape of the expected artifacts follows the approximate contour of the device and extended radially up to 0.8 cm from the implant. When the device is imaged with the sequences prescribed in the standard, the stent lumen cannot be visualized.

#### DRUG INFORMATION

##### Action Mechanism

Rapamycin can act as a cytostatic agent, blocking cell proliferation early in the cell cycle (G1 phase) and returning the cells to their resting state. Rapamycin first binds a cyclophilin FKBP12, and this complex binds and inhibits the function of mTOR (mammalian target of rapamycin), a serine threonine kinase possessing roles such as cellular growth, proliferation, and survival. mTOR acts as a central controller, regulating translation initiation through the eukaryotic initiation factor 4E, and ribosomal p70 S6 kinase pathway.

##### Metabolism

Rapamycin is extensively metabolized by the hepatic CYP3A4 system and is also a substrate in the p-glycoprotein pump of the intestinal wall. The clearance of rapamycin is affected by both of these pathways and displays wide inter-patient variability. The majority of the seven metabolites are formed via O-demethylation or hydroxylation. The metabolites account for less than 10% of the immunosuppressant activity of rapamycin. These metabolites are excreted in bile and feces. The half-life of rapamycin ranges from 57-62 hours, making once daily dosing feasible. Pharmacologic studies have revealed that plasma levels above 15 ng/mL are easily achievable and were well tolerated in patients treated with rapamycin. However, since rapamycin is very lipid soluble, almost no drug is released into the bloodstream during stent placement at the lesion site, and after stent implantation, the diffusion gradient favors elution into tissue again limiting the quantities of circulating free rapamycin.

##### Drug Interactions Following the Oral Administration of Rapamycin

Drug interaction studies have not been conducted with the CORACTO™ Rapamycin-Eluting Coronary Stent Delivery System. Rapamycin is extensively metabolized by the cytochrome P450 3A4 (CYP3A4) in the gut wall and liver and undergoes efflux from enterocytes of the small intestine by P-glycoprotein (P-gp). Therefore, the absorption and subsequent elimination of systemically absorbed rapamycin may be influenced by drugs that affect these proteins. Inhibitors of CYP3A4 and P-gp may increase rapamycin levels, while inducers of CYP3A4 and P-gp may decrease rapamycin levels. The pharmacokinetic interaction between orally-administered rapamycin and concomitantly-administered drugs is discussed below. Drug interaction studies have not been conducted with drugs other than those described below.

##### • Ketoconazole

Multiple-dose ketoconazole administration significantly affected the rate and extent of absorption and Rapamycin exposure after administration of a rapamycin oral formulation, as reflected in the increases in rapamycin

$C_{max}$ ,  $t_{max}$ , and AUC of 4.3-fold, 38%, and 10.9-fold, respectively. However, the terminal  $t_{1/2}$  of rapamycin was not changed. Single-dose rapamycin did not affect steady-state 12-hour plasma ketoconazole concentrations. It is recommended that rapamycin oral solution and oral tablets should not be administered with ketoconazole.

##### • Rifampin

Pre-treatment of 14 healthy volunteers with multiple doses of rifampin, 600 mg daily for 14 days, followed by a single 20 mg dose of rapamycin, greatly increased rapamycin oral-dose clearance by 5.5-fold (range = 2.8 to 10), which represents a mean decrease in AUC and  $C_{max}$  of about 82% and 71%, respectively. In patients where rifampin is indicated, alternative therapeutic agents with lower enzyme induction potential should be considered.

##### • Diltiazem

The simultaneous oral administration of 10 mg of a rapamycin oral solution and 120 mg of diltiazem to 18 healthy volunteers significantly affected the bioavailability of rapamycin. Rapamycin  $C_{max}$ ,  $t_{max}$ , and AUC were increased 1.4-, 1.3-, and 1.6-fold, respectively. Rapamycin did not affect the pharmacokinetics of either diltiazem or its metabolites desacetyldiltiazem and desmethyl diltiazem.

##### • Cyclosporine

Single-dose pharmacokinetic interactions between cyclosporine and rapamycin were investigated for two rapamycin oral formulations in studies using 24 healthy volunteers. Compared to the results obtained when oral rapamycin was administered alone, the oral administration of 10 mg rapamycin 4 hours after a single dose of 300 mg cyclosporine soft gelatin capsules increased mean rapamycin AUC by 33% to 80% and increased mean rapamycin  $C_{max}$  by 33% to 58%, depending on the sirolimus formulation. The half-life of rapamycin was not significantly affected. The cyclosporine mean AUC and mean  $C_{max}$  were not significantly affected.

##### • Drugs which may be co-administered without dose adjustment

Clinically significant pharmacokinetic drug-drug interactions were not observed in studies of the drugs listed below in conjunction with orally-administered rapamycin. Rapamycin and these drugs may be co-administered without dose adjustment.

- Acyclovir
- Digoxin
- Glyburide
- Nifedipine
- Norgestrel/ethinyl estradiol
- Prednisolone
- Sulfamethoxazole/trimethoprim

##### • Other drug interactions

Drugs that may increase rapamycin blood concentrations include:

- Calcium channel blockers: nicardipine, verapamil
- Antifungal agents: clotrimazole, fluconazole, itraconazole
- Macrolide antibiotics: clarithromycin, erythromycin, troleanomycin
- Gastrointestinal prokinetic agents: cisapride, metoclopramide
- Other drugs: bromocriptine, cimetidine, danazol, HIV-protease inhibitors (e.g., ritonavir, indinavir)

Drugs that may decrease rapamycin levels include:

- Anticonvulsants: carbamazepine, phenobarbital, phenytoin
- Antibiotics: rifabutin, rifampin

These lists are not exhaustive.

Care should be exercised when drugs or other substances that are metabolized by CYP3A4 are administered concomitantly with the implantation of CORACTO™ Stents.

##### • Grapefruit juice:

Grapefruit juice reduces the CYP3A4-mediated metabolism of rapamycin.

##### • Herbal Preparations:

St. John's Wort (*Hypericum perforatum*) induces CYP3A4 and P-glycoprotein. As rapamycin is a substrate for both cytochrome

CYP3A4 and P-glycoprotein, there is the possibility that the use of St. John's Wort in patients receiving CORACTO™ Stents may result in reduced rapamycin levels.

##### • Vaccination

Immunosuppressants may affect the response to vaccination. Vaccinations may therefore be less effective for some time after receiving a CORACTO™ Stent. The use of live vaccines should be avoided; live vaccines may include, but are not limited to, those for measles, mumps, rubella, oral polio, BCG, yellow fever, varicella, and TY21a typhoid.

##### • Drug-laboratory test interactions

There have been no studies on the interactions of rapamycin in commonly employed clinical laboratory tests.

##### Mutagenesis, Carcinogenicity and Reproductive Toxicology

The genotoxicity, carcinogenicity and reproductive toxicity of CORACTO™ Stents have not been evaluated. However, the genotoxicity, carcinogenicity, and reproductive toxicity of rapamycin have been investigated in bacterial and mammalian cells in vitro and in laboratory animals in vivo.

Rapamycin was not genotoxic in the in vitro bacterial reverse mutation assay, the Chinese hamster ovary cell chromosomal aberration assay, the mouse lymphoma cell forward mutation assay, or the in vivo mouse micronucleus assay.

Potential Adverse Events Related to Rapamycin (Following Oral Administration):

- Abnormal liver function tests
- Anemia
- Arthralgias
- Diarrhea
- Hypercholesterolemia
- Hypersensitivity, including anaphylactic/anaphylactoid-type reactions
- Hypertriglyceridemia
- Hypokalemia
- Infections
- Interstitial lung disease
- Leukopenia
- Lymphoma and other malignancies
- Thrombocytopenia

##### POTENTIAL ADVERSE EVENTS

- Additional intervention (Percutaneous, Surgical Intervention (Emergency CABG))
- Air Embolism
- Allergic or drug reactions to antiplatelet agents /anticoagulation agents/ contrast medium/stent material
- Arrhythmia
- Cardiac Tamponade
- Damage to healthy intima
- Death
- Dissection
- Hemorrhage/Bleeding
- Hematoma
- Infection
- Myocardial Infarction
- Myocardial Ischemia
- Perforation
- Pseudoaneurysm
- Renal Failure
- Restenosis (Requiring and/or not requiring re-intervention)
- Stent Embolism
- Stroke
- Thrombosis
- Total Occlusion of the coronary artery

- Vessel Spasm

#### INSTRUCTIONS FOR USE

Materials Required (Note: These materials are not included in the package of the product)

MATERIAL
Appropriate guiding catheter
10-20 cc syringes
Sterile Heparinized Saline
0.014" guidewire
Introducer sheath
Rotating hemostatic valve
Contrast diluted 1:1 with normal saline
Inflation device
Three-way stopcock
Appropriate anticoagulation and antiplatelet medication

#### A. PREPARATION

1. Inspect the package of the product against the damage of the packaging and sterility barrier.

##### Notes:

- Do not use the device in case inner package sterility is compromised, package is damaged or device labeling is incomplete.
- Before starting the procedure, prepare the inflation device according to the Instructions for Use.

2. Carefully remove the delivery system from its protective tubing.
3. Ensure that the diameter and length of the stent is appropriate according to the treated vessel.

**Note:** Correct stent sizing is important for successful stenting. In general stent size should be chosen to match the diameter of the reference vessel and to correspond to the length of the lesion. Slight stent oversizing is preferable to undersize.

4. Gently remove the product from the packaging and carefully inspect for kink, bend and any type of damages. Avoid contacting any foreign surfaces before using in the patient to maintain the sterility of the product.
5. With the tip down, orient stent delivery system vertically. Examine the product and verify that the stent is centered between radiopaque gold markers.

##### Notes:

- Any damage to delivery system may decrease the desired performance characteristics.
  - Do not use if any defects are noted.
  - Do not touch or handle the stent in order to prevent stent damage or stent dislodgement.
6. Prepare inflation device with diluted contrast medium.
  7. Fill the 20 cc syringe with saline solution. Carefully flush carefully the guidewire lumen in routine manner.
  8. Attach inflation device to a stopcock, and then attach the stopcock to the inflation port. Make sure that the components are connected appropriately. Prepare the inflation device in accordance with the manufacturer's instructions.
  9. After delivering the stent in the appropriate position at the target lesion open the stopcock to the delivery system and pull negative until de-airing the inflation lumen of the delivery system and then release to neutral position.

##### Notes:

- Do not use air or any gaseous medium to inflate the balloon.
  - Do not pull negative pressure on the stent delivery system before positioning the stent in the target lesion.
10. This procedure (step 9) can be repeated 2 to 3 times to completely de-air the stent delivery system.

##### Notes:

- The same procedure can be done with a syringe instead of the

inflation device.

- Do not apply negative pressure on the system with the inflation device while introducing the delivery system into vessel. This may cause stent dislodgement.

11. Ensure catheter sheath introducer/guiding catheter stability.

#### B. DELIVERY AND DEPLOYMENT PROCEDURE

1. Prepare vascular access site according to standard PTCA practice.
2. It is highly recommended to pre-dilate the lesion with a PTCA catheter.
3. Advance the stent delivery system over the guidewire (max. diam. 0.014") through the hemostatic valve. The hemostatic valve should be completely open while the product is introduced.

**Note:** If you encounter any resistance, do not force passage, as this may cause stent damage or stent dislodgement from the balloon. Ensure that the guiding catheter is in the correct position and that its tip is stationary.

4. Advance the stent delivery system under direct fluoroscopic visualization.
5. Position the stent across the lesion using the proximal and distal radiopaque markers on the balloon as the reference point.
6. To expand the stent gently inflate the balloon to the targeted vessel diameter. Follow the product labeling for the proper inflation pressure. The fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum stent diameter compared to the vessel. Optimal expansion and proper sizing requires that the stent is full in contact with the arterial wall.

##### Notes:

- Expansion of the stent should not be undertaken if the stent positioning is not optimal. Do not exceed Rated Burst Pressure as indicated on the label. Stent overexpansion and under expansion may lead to several complications
  - Do not expand the stent beyond the compliance chart provided on the label.
  - Should any resistance be felt at any time during either lesion access or removal of delivery system post-stent implantation, the entire system should be removed as a single unit.
7. Tighten the rotating hemostatic valve.

#### C. DEPLOYMENT AND REMOVAL PROCEDURE

1. Deflate the balloon by pulling negative on inflation device. Give adequate time for full-balloon deflation. Deflation of the balloon should be confirmed by absence of contrast within the balloon.
2. Fully open rotating hemostatic valve.
3. Very slowly withdraw the balloon from the stent maintaining negative pressure to gently dislodge the balloon from the stent.
4. Once the balloon has been completely deflated, carefully withdraw the delivery system from the guiding catheter under fluoroscopy while maintaining the guidewire positioned.
5. Tighten the rotating hemostatic valve.

**Note:** If any unusual resistance is encountered when withdrawing the delivery system, the stent delivery system and guiding catheter should be removed as a single unit. This must be done under direct visualization with fluoroscopy.

6. Repeat angiographic imaging to assess the stented area. If the stent is under-dilated, carefully position the original delivery system or another catheter with appropriate diameter and repeat the previous steps to expand the stent to the required diameter.

Stent retrieval methods (use of additional wires, snares and / or forceps) may result in additional trauma to the coronary vessel and / or the vascular access site. Complications may include bleeding, hematoma or pseudoaneurysm.

#### TERMS OF WARRANTY

Alvimedica warrants that reasonable care has been used in the design and manufacture of this instrument. Due to biological variations among individuals, no product is 100 % effective in every case. Therefore, and since Alvimedica has no control over conditions under which the product used after sales, the selection of patients and the methods of application; Alvimedica disclaims all warranties with respect to the product, expressed or implied including, without limitation, those of merchantability and/or fitness for a particular purpose. Alvimedica is not directly or indirectly responsible for any injury or damage to or loss of any person resulting from the misuse of the product, nor is Alvimedica responsible or liable, directly or indirectly, for any injury, damage or loss that may result from re-use or re-sterilization.

#### INTELLECTUAL PROPERTY PROTECTION

All trademarks, logos, inventions, know-how, technology, proprietary information and other intellectual property rights of this product solely belongs to Alvimedica. Alvimedica reserves all the rights to demand, sue, and/or seek relief before judicial authorities in case of an infringement regarding any unauthorized use or misuse of intellectual property or reverse engineering related the product.

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